


Magnetic Disc Instructions

This is a flexible magnetic rubber sheet, which when mounted to a platen, is then able to hold all magnetic versions of polishing pads and grinding discs. If used with a thin metal plate (KMS), it can be used with self adhesive grinding papers. This is a magnet based system.

1. Peel off the brown barrier. This will reveal a sticky layer



2. Stick the black magnetic disc to the platen



Kemet International Limited

Parkwood Trading Estate, Maidstone, Kent, ME15 9NJ, UK
Tel: +44(0) 1622 755287; Fax: +44(0) 1622 670915
Email: sales@kemet.co.uk; Web: www.kemet.co.uk

3. The magnetic platen can be used with Kemet metallic polishing pads, or magnetic consumables from other suppliers



Kemet International Limited

Parkwood Trading Estate, Maidstone, Kent, ME15 9NJ, UK
Tel: +44(0) 1622 755287; Fax: +44(0) 1622 670915
Email: sales@kemet.co.uk; Web: www.kemet.co.uk